
AH Serisi Otomatik Homojenizatör

İçerik

01

Homojenizasyon

02

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03

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04

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05

Cihaz Konfigürasyonu

01

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1.1 Introduction of AH-50 and AH-40

AH-50
Automated Homogenizer

AH-40
Automated Homogenizer

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Explain:Uniform when the dispersion tool down to the most low stop time.

1

2

3

4

RayKol Ready

+ New Open Save Save as

Home

Method

Configure

Maintenance

Help

Name	Speed	Volume	Time	Channel Sample	Other
Filling liquid					
Filling liquid					
Filling liquid					
Homogenize					
Water					
Ultrasonic					
Solvent					
End					

Pre filling

Pre filling setting

Note:
The water rinsing setting range is 0-1, pre-setting is "0". When "0" the cutter will be without up and down moving and rotating, and "1" the cutter will run up and down, and rotate. Note: The water rinsing setting range is 0-1, pre-setting is "0". When "0" the cutter will be without up and down moving and rotating, and "1" the cutter will run up and down, and rotate.

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
GC, GCMS, LCMS

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Gıdalardaki herbisitler:


- Determination of cyclohexanone herbicides in food by LC-MS-MS
- Determination of multiple ether herbicides residues in foods by GC-MS
- Determination of thiocarbamate herbicides in foods by LC-MS-MS
- Determination of amitrole herbicides in foods by GC-MS
- Determination of Dinoseb Residues in Food by LC-MS

Gıdalardaki pestisitler:

- Determination of thiamethoxam and its metabolite thiamethoxam residues in food by LC-MS-MS
- Determination of diflubenzuron residues in food by LC-MS-MS
- Determination of 9 carbamate pesticides and their metabolite residues in plant-derived foods by LC-Post-column derivatization
- Determination of 90 organophosphorus pesticides and metabolites residues in plant-derived foods
- Determination of abamectin residues in food by LC-MS-MS

Çay yaprağındaki pestisitler:

- Determination of residues of 448 pesticides and related chemicals in tea by LC-MS

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